CONSUMER ADVANTAGES

With Platinum HD Mobile, wearers take advantage of these benefits:

REDUCED VISUAL FATIGUE

Easy and agile switch from the

DYNAMIC VISION

near to intermediate zone

IMPROVED FOCUS Better ability to read text

on digital displays

NATURAL POSITION Ergonomically comfortable

More relaxed eyes

FASTER ADAPTATION Almost immediate adaptation **FULLY PERSONALIZED**

Unique lenses for each patient

PERSONALIZATION PARAMETERS

To provide 100% personalization, it is essential to include all personalization parameters unique to each wearer's prescription information. For orders that do not include personalization parameter data (i.e. pantoscopic angle, wrap angle, monocular PD, back vertex distance, and frame dimensions), the lens will be optimized using default values.

MINIMUM FITTING HEIGHTS

14 mm, 15 mm, 16 mm, 17 mm, 18 mm

MOUNTING INSTRUCTIONS

For proper positioning of the lens in the frame, the invisible engravings must be considered.



800.634.9461 • mjoptical.com 8838 Washington Circle, Omaha, NE 68127


PLATINUM HD mobile



A fully personalized progressive lens ideal for **electronic device** users









OPTIMUM RESOLUTION THANKS TO DIGITAL RAY-PATH® TECHNOLOGY

Digital Ray-Path® is an innovative calculation technology that uses a design engine to compensate the lens with a simulation of the binocular eye-lens system. Every unique lens is individually calculated guaranteeing an adapted solution for any prescription and base curve.



SMART ADD TECHNOLOGY ALLOWS THE EYES TO BE RELAXED WHEN VIEWING SCREEN DISPLAYS

Platinum HD Mobile incorporates Smart Add, the technology that improves the patient's visual experience when using digital devices. This technology assists with changes in focus at different working distances in a more agile and accurate way when the patient is working or reading simultaneously with different screens.

TARGET & POSITIONING

• Premium personalized progressive lens for electronic device users.

• Ideal for progressive lens wearers


ages 40 and over, both experts

SHORT CORRIDOR ALLOWS QUICK TRANSITION TO NEAR VISION

comfortable transition to the reading zone.

Our constant use of smartphones and tablets has led to a greater need for switching between near and distance vision. Keeping the same corridor length, Platinum HD Mobile's unique progression profile offers a more comfortable transition from the distance zone to the reading zone.

PROGRESSION PROFILE



SMOOTH PROGRESSION PROFILE FOR GREATER COMFORT

The transition between distance and near vision is very smooth. Thanks to its unique layout, the Platinum HD Mobile lens offers a wide and comfortable reading zone.















and novices.





Variable available Inset

Digital

personalized Ray-Path®

Smart Add

Multiple corridor

Short corridor available